


Glencoe Science

LAB MANUAL

Probeware for Biology



science.glencoe.com



Glencoe

New York, New York Columbus, Ohio Chicago, Illinois Peoria, Illinois Woodland Hills, California

Credits

Each *Probeware Lab Manual for Biology* activity was reviewed by Vernier Software & Technology.



The McGraw-Hill Companies

Copyright © by The McGraw-Hill Companies, Inc. All rights reserved. Except as permitted under the United States Copyright Act, no part of this publication may be reproduced or distributed in any form or by any means, or stored in a database retrieval system, without prior written permission of the publisher.

The terms CBL 2, TI-GRAPH LINK, TI Connect and TI Interactive! are either registered trademarks of, trademarks of, or copyrighted by Texas Instruments, Inc. Vernier LabPro is a registered trademark of and Graphical Analysis is copyrighted by Vernier Software & Technology. PASPORT is a registered trademark of PASCO Scientific. Macintosh is a registered trademark of Apple Computer, Inc. Windows is a registered trademark of Microsoft Corporation in the United States and/or other countries.

Send all inquiries to:
Glencoe/McGraw-Hill
8787 Orion Place
Columbus, OH 43240-4027

ISBN 0-07-860226-2

Printed in the United States of America.

1 2 3 4 5 6 7 8 9 10 009 08 07 06 05 04 03

Contents

To the Student	iv
Getting Started with Probeware	v
Laboratory and Safety Guidelines	vii
Safety Symbols	viii
Lab 1 Is oxygen cycled in the environment?	1
Lab 2 An Environmental Limiting Factor	5
Lab 3 Testing Water Quality	9
Lab 4 How well does yeast ferment different sugars?	13
Lab 5 How can pH be used to compare rates of photosynthesis?	17
Lab 6 Effect of Environmental Temperature on the Metabolic Rates of Animals	21
Lab 7 What is the effect of exercise on body temperature?	25
Lab 8 Measuring Response Time	29
Lab 9 Breathing and Heart Rate	33
Lab 10 What is the effect of exercise on heart rate?	37
Appendix	41

To the Student

The activities in this book are designed to help you study science using probeware technology. A probeware lab is different from other labs because it uses a probe or sensor to collect data, a data collection unit to interpret and store the data, and a graphing calculator or computer to analyze the data. These components are connected with a software program called DataMate that makes them work together in an easy-to-use, handheld, system. These labs are designed specifically for the TI-73 or TI-83 Plus graphing calculators and a CBL 2™ (produced by Texas Instruments, Inc.) or Vernier LabPro® (produced by Vernier Software & Technology) data collection unit.

The activities in this book will help you improve your ability to recognize and use equipment properly and to analyze data. To help you get started, the next few pages will provide you with:

- information about **getting started with probeware**
- a list of **laboratory and safety guidelines**
- a reference page of **safety symbols**

Each lab activity in this manual includes the following sections:


- **Introduction**—The introductory paragraphs give you background information needed to understand the activity.
- **Objectives**—The list of objectives is a guide to what will be done in the activity and what will be expected of you.
- **Materials**—The materials section lists the supplies you will need to complete the activity.
- **Procedure**—The procedure gives you step-by-step instructions for carrying out the activity. Many steps have safety precautions. Be sure to read these statements and obey them for your own and your classmates' protection. Unless told to do otherwise, you are expected to complete all parts of each assigned activity. Important information needed for the procedure but that is not an actual procedural step is also found in this section.
- **Hypothesis**—You will develop a hypothesis statement to express your expectations of the results and as an answer to the problem statement.
- **Data and Observations**—This section includes tables and space to record data and observations.
- **Analysis**—In this section, you draw conclusions about the results of the activity just completed. Rereading the introduction before answering the questions is most helpful at this time.
- **Checking Your Hypothesis**—You will determine whether your data supports your hypothesis.
- **Further Investigations/Further Explorations**—This section gives ideas for further activities that you may do on your own. They may be either laboratory or library research.

Getting Started with Probeware

The following instructions will guide you through the setup process for the data collection unit and the graphing calculator. The activities are compatible with either the CBL 2 or the LabPro unit. Each activity was written for use with TI-73 or TI-83 Plus graphing calculators. These activities can be adapted for use with other graphing calculators or other data collection units, if desired.

Connecting a Graphing Calculator to the CBL 2 or LabPro Unit

1. Insert batteries into the CBL 2 or LabPro unit and graphing calculator.
2. The cradle is an optional accessory that conveniently connects the two units. Slide the back of the cradle onto the front of the CBL 2 or LabPro unit until it clicks into place.
3. Insert the upper end of the calculator into the cradle and press down on the lower end until it locks into place.
4. Connect the CBL 2 or LabPro unit to the graphing calculator using the unit-to-unit link cable. Plug the cable into the I/O port at the end of the CBL 2 or LabPro unit and the other end into the I/O port at the end of the calculator. Make sure that the unit-to-unit link cable is securely in place.



Resetting the Calculator Memory

It is recommended that the memory of the calculator be cleared before the DataMate data collection program is transferred.

1. Press $\boxed{2\text{nd}}$ [MEM].
2. Select **Reset**.
3. Select **ALL RAM . . .**
4. Select **Reset**. The calculator screen will display **RAM cleared**.

Transferring DataMate to the Calculator

The DataMate program is stored on the CBL 2 or LabPro unit and is transferred to the graphing calculator for use. Once DataMate is transferred to the graphing calculator, it will remain there until the calculator memory is reset using the instructions above.

1. For the TI-73, press $\boxed{\text{APPS}}$. Select **Link . . .**
For the TI-83 Plus, press $\boxed{2\text{nd}}$ [LINK].
2. Use the right arrow to highlight **RECEIVE**. Press $\boxed{\text{ENTER}}$.
3. The screen will display **Waiting . . .** Press the large **TRANSFER** key found on the upper left-hand side of the CBL 2 or LabPro unit. When the transfer is complete, the screen will display the transferred programs followed by the word **Done**.
4. Press $\boxed{2\text{nd}}$ [QUIT].

Starting DataMate

When you are ready to collect data, use the following instructions to start DataMate.

For the TI-73:

1. Press $\boxed{\text{PRGM}}$.
2. Select **DataMate**.
3. Press $\boxed{\text{ENTER}}$.

For the TI-83 Plus:

1. Press $\boxed{\text{APPS}}$.
2. Select **DataMate**.

Setting up Probes Manually

The CBL 2 and LabPro unit should recognize the probe attached automatically. If this does not happen, follow these instructions.

1. Select **SETUP** from the DataMate main screen.
2. Press $\boxed{\text{ENTER}}$ to select channel 1, or press $\boxed{\blacktriangle}$ to select the channel where the probe is inserted and then press $\boxed{\text{ENTER}}$.
3. Select the correct sensor number from the SELECT SENSOR menu.
4. If requested, select the type of probe used.
5. Select **OK** to return to the DataMate main screen.

Laboratory and Safety Guidelines

Emergencies

- Inform the teacher immediately of *any* mishap—fire, injury, glassware breakage, chemical spills, and so forth.
- Know the location of the fire extinguisher, safety shower, eyewash, fire blanket, and first-aid kit. Know how to use this equipment.
- If chemicals come into contact with your eyes or skin, flush with large quantities of water and notify your teacher immediately.

Preventing Accidents

- Do NOT wear clothing that is loose enough to catch on anything. Do NOT wear sandals or open-toed shoes. Remove loose jewelry—chains or bracelets—while doing lab work.
- Wear protective safety gloves, goggles, and aprons as instructed.
- Always wear safety goggles (not glasses) in the laboratory.
- Wear goggles throughout the entire activity, cleanup, and handwashing.
- Keep your hands away from your face while working in the laboratory.
- Remove synthetic fingernails before working in the lab (these are highly flammable).
- Do NOT use hair spray, mousse, or other flammable hair products just before or during laboratory work where an open flame is used (they can ignite easily).
- Tie back long hair and loose clothing to keep them away from flames and equipment.
- Eating, drinking, chewing gum, applying makeup, and smoking are prohibited in the laboratory.
- Do NOT inhale vapors or taste, touch, or smell any chemical or substance unless instructed to do so by your teacher.

Working in the Laboratory

- Study all instructions before you begin a laboratory or field activity. Ask questions if you do not understand any part of the activity.
- Work ONLY on activities assigned by your teacher. NEVER work alone in the laboratory.
- Do NOT substitute other chemicals/substances for those listed in your activity.
- Do NOT begin any activity until directed to do so by your teacher.
- Do NOT handle any equipment without specific permission.
- Remain in your own work area unless given permission by your teacher to leave it.
- Do NOT point heated containers—test tubes, flasks, and so forth—at yourself or anyone else.
- Do NOT take any materials or chemicals out of the classroom.
- Stay out of storage areas unless you are instructed to be there and are supervised by your teacher.


Laboratory Cleanup


- Keep work, lab, and balance areas clean, limiting the amount of easily ignitable materials.
- Turn off all burners, water faucets, probeware, and calculators before leaving the lab.
- Carefully dispose of waste materials as instructed by your teacher.
- With your goggles on, wash your hands thoroughly with soap and warm water after each activity.


Safety Symbols


SAFETY SYMBOLS

	HAZARD	EXAMPLES	PRECAUTION	REMEDY
DISPOSAL 	Special disposal procedures need to be followed.	certain chemicals, living organisms	Do not dispose of these materials in the sink or trash can.	Dispose of wastes as directed by your teacher.
BIOLOGICAL 	Organisms or other biological materials that might be harmful to humans	bacteria, fungi, blood, unpreserved tissues, plant materials	Avoid skin contact with these materials. Wear mask or gloves.	Notify your teacher if you suspect contact with material. Wash hands thoroughly.
EXTREME TEMPERATURE 	Objects that can burn skin by being too cold or too hot	boiling liquids, hot plates, dry ice, liquid nitrogen	Use proper protection when handling.	Go to your teacher for first aid.
SHARP OBJECT 	Use of tools or glassware that can easily puncture or slice skin	razor blades, pins, scalpels, pointed tools, dissecting probes, broken glass	Practice common-sense behavior and follow guidelines for use of the tool.	Go to your teacher for first aid.
FUME 	Possible danger to respiratory tract from fumes	ammonia, acetone, nail polish remover, heated sulfur, moth balls	Make sure there is good ventilation. Never smell fumes directly. Wear a mask.	Leave foul area and notify your teacher immediately.
ELECTRICAL 	Possible danger from electrical shock or burn	improper grounding, liquid spills, short circuits, exposed wires	Double-check setup with teacher. Check condition of wires and apparatus.	Do not attempt to fix electrical problems. Notify your teacher immediately.
IRRITANT 	Substances that can irritate the skin or mucous membranes of the respiratory tract	pollen, moth balls, steel wool, fiberglass, potassium permanganate	Wear dust mask and gloves. Practice extra care when handling these materials.	Go to your teacher for first aid.
CHEMICAL 	Chemicals that can react with and destroy tissue and other materials	bleaches such as hydrogen peroxide; acids such as sulfuric acid, hydrochloric acid; bases such as ammonia, sodium hydroxide	Wear goggles, gloves, and an apron.	Immediately flush the affected area with water and notify your teacher.
TOXIC 	Substance may be poisonous if touched, inhaled, or swallowed.	mercury, many metal compounds, iodine, poinsettia plant parts	Follow your teacher's instructions.	Always wash hands thoroughly after use. Go to your teacher for first aid.
OPEN FLAME 	Open flame may ignite flammable chemicals, loose clothing, or hair.	alcohol, kerosene, potassium permanganate, hair, clothing	Tie back hair. Avoid wearing loose clothing. Avoid open flames when using flammable chemicals. Be aware of locations of fire safety equipment.	Notify your teacher immediately. Use fire safety equipment if applicable.

 **Eye Safety**
Proper eye protection should be worn at all times by anyone performing or observing science activities.

 **Clothing Protection**
This symbol appears when substances could stain or burn clothing.

 **Animal Safety**
This symbol appears when safety of animals and students must be ensured.

 **Radioactivity**
This symbol appears when radioactive materials are used.

Lab 1

Is oxygen cycled in the environment?

Probeware Activity

Plants and animals interact with each other in many ways. One of those ways is by cycling carbon. For example, animals release the carbon dioxide that is produced when their cells break down food. Plants take in the carbon dioxide and use it to make food during photosynthesis. One of the by-products of photosynthesis is oxygen. What happens to the oxygen that is produced by plants? Does it cycle between plants and animals? In this lab, you will design an experiment to find out, using a probe that measures the concentration of oxygen dissolved in water.

Problem

How could you find out whether oxygen is cycled between plants and animals in the environment? What would happen if you placed an aquatic plant and/or aquatic animal in a closed aquatic system, one in which gases could not enter or leave?

Hypothesis

Hypothesize whether oxygen is cycled between plants and animals in the environment. Write your hypothesis below.

Objectives

- Hypothesize whether oxygen is cycled between plants and animals in the environment.
- Design a closed aquatic system that does not allow gases to enter or leave.
- Interpret data to determine whether oxygen is cycled between plants and animals.

Possible Materials

- | | | |
|---|---|--|
| <input type="checkbox"/> LabPro or CBL 2 interface | <input type="checkbox"/> aquatic animal, such as a snail | <input type="checkbox"/> 100-mL graduated cylinder |
| <input type="checkbox"/> TI graphing calculator | <input type="checkbox"/> aquatic plant, such as <i>Elodea</i> | <input type="checkbox"/> 10-mL graduated cylinder |
| <input type="checkbox"/> LabPro or CBL 2 AC adapter (optional) | <input type="checkbox"/> tap water (allowed to stand for one day) | <input type="checkbox"/> 250-mL beakers |
| <input type="checkbox"/> link cable | <input type="checkbox"/> artificial light source | <input type="checkbox"/> distilled water |
| <input type="checkbox"/> Vernier dissolved oxygen probe | <input type="checkbox"/> D.O. electrode filling solution | <input type="checkbox"/> wax marking pencil |
| <input type="checkbox"/> small, clear containers with caps, to be used for the closed systems | <input type="checkbox"/> Beral pipette | <input type="checkbox"/> lab wipes |
| | <input type="checkbox"/> metric ruler | <input type="checkbox"/> laboratory apron |
| | | <input type="checkbox"/> goggles |

Lab

1

Is oxygen cycled in the environment?,

continued

Probeware Activity

Plan the Experiment

1. Devise a procedure you can use to determine whether oxygen is cycled between plants and animals in an aquatic environment. The procedure should use a probeware system that includes a dissolved oxygen probe, LabPro or CBL 2 interface, TI graphing calculator, and link cable.
2. Decide on the type of closed system you will design and the number of closed systems you will need to observe in order to make valid conclusions. Decide on the variable that you will manipulate and those that you will keep constant.
3. Think about how long you will wait before taking measurements of the dissolved oxygen levels in the closed systems. Before you collect data, the probeware system including the dissolved oxygen probe must be set up, warmed up, and calibrated. Allow time in your experimental plan for these steps.
4. Write your procedure on another sheet of paper or in your notebook. It should include all the materials you will use.

CAUTION: *To avoid harming the plants and animals, do not place the closed systems too close to an artificial light source. Heat from the light source may raise the water temperature above a safe level.*

Check the Plan

1. Be sure to include a control group for comparison purposes in your experiment. Do you have different combinations of organisms in the experimental groups?
2. Make sure the teacher has approved your experimental plan before you proceed further.

Carry Out Your Experiment

1. Follow the steps in your plan.
2. Connect the TI graphing calculator to the LabPro or CBL 2 interface using the link cable. Connect the dissolved oxygen probe

into Channel 1 of the interface. If the dissolved oxygen probe needs to be warmed up, proceed to Step 3. If the probe has already been warmed up, proceed to Step 7.




Figure 1

3. Set up the dissolved oxygen probe.
 - a. Unscrew the membrane cap (counterclockwise) from the tip of the electrode on the dissolved oxygen probe. Do not touch the membrane at the very tip of the probe.
 - b. Use a Beral pipette to fill the membrane cap with about 1 mL of D.O. electrode filling solution. Carefully thread the membrane cap (clockwise) onto the electrode body. Do not over tighten the cap. Rinse the electrode with distilled water and carefully wipe it dry with a lab wipe.
 - c. Place the dissolved oxygen probe in a 250-mL beaker containing about 75 mL of distilled water.
4. Turn on the calculator and start the DATA-MATE program. Press **CLEAR** to reset the program.
5. The dissolved oxygen probe must be powered and connected to the LabPro or CBL 2 interface to be warmed up.
 - a. If the calculator screen displays CH 1 DO (MG/L), proceed to Step 6. If it does not, continue with this step to manually select the dissolved oxygen probe.

Lab

1

Is oxygen cycled in the environment?,

continued

Probeware Activity

- b. Select SETUP from the main screen.
 - c. Press **ENTER** to select CH 1.
 - d. Select D. OXYGEN (MG/L) from the SELECT SENSOR menu.
 - e. Select OK to return to the main screen.
6. Warm up the dissolved oxygen probe for 10 minutes.
- a. With the probe still in the water, wait 10 minutes while the probe warms up. The probe must stay connected to the interface at all times to keep it warmed up. If disconnected for a period longer than a few minutes, it will be necessary to warm it up again.
 - b. At the end of class, leave the dissolved oxygen probe connected to the interface, with the DATAMATE program running. If this is done, the probe will stay warm and ready for the next class.
7. **NOTE:** *Proceed with this part of the lab only after the dissolved oxygen probe has been warmed up.*
- a. After the desired number of days has passed, open one of the closed systems, being careful not to agitate the water, which could cause a change in the amount of dissolved oxygen.
 - b. Slowly pour the water from the closed system into a 250-mL beaker, and then insert the dissolved oxygen probe into the water. Using a gentle swirling motion, stir the dissolved oxygen probe through the water. Make sure no bubbles are trapped under the tip of the probe. Liquid must be continually moving past the membrane of the electrode. Monitor the dissolved oxygen readings displayed on the calculator screen. Once the reading has stabilized, record the value in Table 1, or you can make your own data table. Rinse the end of the probe with distilled water before testing the water from another closed system.
- c. Place the probe in a beaker of distilled water. Leave the DATAMATE program and calculator running for the next class.
 - d. If you are the last class to use the equipment, exit the DATAMATE program and turn off the calculator. Disconnect the probe from the LabPro or CBL 2. Remove the membrane cap and rinse the inside and outside of the cap with distilled water. Rinse and carefully dry the exposed cathode and anode inner elements of the probe. Reinstall the membrane cap loosely onto the electrode body for storage.
8. After doing the experiment, return the organisms to their original locations or follow the teacher's directions.

Data and Observations

Table 1

Descriptions of Closed Systems	Concentration of Dissolved Oxygen (ppm)

Lab

1

Is oxygen cycled in the environment?,*continued***Probeware Activity****Analysis**

1. Which closed system had the highest concentration of dissolved oxygen? Explain why.

2. Which closed system had the lowest concentration of dissolved oxygen? Explain why.

3. From your data, what can you conclude about whether oxygen is cycled between plants and animals? Explain.

4. Why was it necessary to have a closed system in your experiment?

Checking Your Hypothesis

Was your hypothesis supported by your data? Why or why not?

Further Investigations

1. Repeat the lab using different aquatic plants and/or animals. Compare your results from the two experiments.
2. Modify the lab to compare the dissolved oxygen levels for closed systems that have been placed in low-light environments or in the dark.

Lab 2

An Environmental Limiting Factor

Probeware Activity

Just as canaries warned early coal miners of dangerous levels of gas, so are fish, plants, and other organisms important indicators of the health of water and soil environments. These organisms are affected by the acidity of water and soil. The acidity of a substance is measured on the pH scale. Acids have pH values lower than 7, while bases have pH values above 7. A substance with a pH of 7 is neutral, being neither acidic nor basic. The farther a pH value is from 7, the more acidic or basic the substance is.

In this lab, you will investigate pH as a limiting factor by measuring the pH of a variety of soil and water samples in your area. Fish generally do well in pH ranges of 6.7 to 8.5. When the pH drops below 5 or exceeds 9, most fish have difficulty surviving. Thus, pH is considered a limiting factor—an environmental factor that affects the ability of organisms to survive. The acidity of lake waters is affected by acid rain and by minerals that leach out of alkaline (basic) or acidic soils and drain into lakes. Organisms that live in soil also are sensitive to acidity levels. For example, rhododendrons and azaleas are acid loving and do well in soils of pH 4–5.

Objectives

- Use a pH sensor to measure pH.
- Measure and compare the pH of soil and water samples.
- Evaluate the pH of the samples as an environmental limiting factor.

Materials

- | | | |
|---|--|---|
| <input type="checkbox"/> LabPro or CBL 2 unit | <input type="checkbox"/> plastic bags (3 or more) | <input type="checkbox"/> balance |
| <input type="checkbox"/> TI graphing calculator | <input type="checkbox"/> large spoon | <input type="checkbox"/> jars with lids (8) |
| <input type="checkbox"/> link cable | <input type="checkbox"/> distilled water | <input type="checkbox"/> wax marking pencil |
| <input type="checkbox"/> pH sensor | <input type="checkbox"/> water samples (3 or more) | <input type="checkbox"/> plastic or rubber gloves |
| <input type="checkbox"/> 250-mL beakers (2) | <input type="checkbox"/> paper plates (3) | <input type="checkbox"/> laboratory apron |
| <input type="checkbox"/> soil samples (3 or more) | <input type="checkbox"/> mortar and pestle | <input type="checkbox"/> goggles |

Lab

2

An Environmental Limiting Factor, *continued*

Probeware Activity

Procedure

Part A. Preparing Soil Samples

1. Several days before the lab, use a spoon to collect three or more soil samples, each approximately three square centimeters wide by two centimeters deep. Samples could be taken from a lawn, wooded area, wetland, or home garden. Place the samples in separate plastic bags. While collecting, record in Table 1 the kinds of plants growing in the soil. **CAUTION: Wear protective gloves while collecting and handling soil samples.**
2. Break up the samples and let them dry on separate paper plates for several days.
3. Grind the dry soil for each sample with a mortar and pestle. Remove any plant material or rock.
4. For each sample, place about 30 grams of soil and 60 grams of distilled water into a jar. Close the jar, label the source of the soil, and shake the jar vigorously. Allow the water to stand overnight.

Part B. Preparing Water Samples

1. Before the lab, collect water samples in jars with lids. These should include samples of distilled water, ocean or saltwater, aquarium water, precipitation (rain or snow), drinking water, and stream, pond, or lake water. Label each jar by source. Then record the source of each sample in Table 2. **CAUTION: Wear protective gloves while collecting and handling water samples.**
2. Store samples in a cool place. Allow snow or ice to melt at room temperature.

Part C. Preparing the pH Sensor

1. Plug the pH sensor into Channel 1 of the LabPro or CBL 2 interface. Using the link cable, connect the TI graphing calculator to the interface. Push the link cable securely into each jack.

2. Turn on the graphing calculator. Start the DATAMATE program and go to the MAIN MENU. Press **CLEAR** to reset the program. If the DATAMATE program is not loaded, transfer the program from the memory of the LabPro or CBL 2 interface to the TI graphing calculator.
3. Set up the calculator and interface for a pH Sensor.
 - a. If the calculator displays PH in CH 1, proceed directly to Part D. If it does not, continue with this step to set up your sensor manually.
 - b. Select SETUP from the main screen.
 - c. Press **ENTER** to select CH 1.
 - d. Select PH from the SELECT SENSOR menu.
 - e. Select OK to return to the main screen. Readings from the pH sensor will be displayed on the main screen.




Figure 1

Lab

2

An Environmental Limiting Factor, *continued*

Probeware Activity

Part D. Testing the Samples

1. Before each use of the pH sensor, hold it over a beaker and gently rinse the tip with distilled water. **CAUTION: The tip of the pH sensor can be broken easily.** Do not let the sensor dry out. When not in use, immerse the sensor in a beaker of tap water.
2. Pour a soil solution or water sample into a clean, dry beaker.
3. Carefully place the pH sensor in the beaker as shown in Figure 1. Allow the pH value reading to stabilize for 15–20 seconds. Record the value in Table 1 for a soil sample or in Table 2 for a water sample.
4. Dispose of the sample and rinse and dry the beaker.
5. Repeat Part D, steps 1–4, for your other samples.
6. At the conclusion of the lab, wash your hands thoroughly with soap and water.
7. Share your data with other class members.

Data and Observations

Table 1

Source of Soil Sample	Plants That Grow in This Soil	pH

Table 2

Source of Water Sample	pH

Analysis

1. Compare the pH values of the soil and water samples. Which sample showed the highest pH value? The lowest pH value?

2. Did soils of differing pH levels support different types of plants? Explain.

3. Compare the pH values for the precipitation samples and for the samples from bodies of water. Summarize your observations.

4. How does the pH value of your drinking water compare with that of natural bodies of water and precipitation?

5. Which of the water samples you tested would provide the greatest limiting factor for fish?

Further Investigations

1. Contact a nursery or do research to determine how the soil pH in your area acts as a limiting factor for plants.
2. GLOBE (Global Learning and Observations to Benefit the Environment) is a worldwide network of students, teachers, and scientists working together to study and understand the global environment. Join the GLOBE schools to observe the archival data or atmospheric and hydrologic data. GLOBE students make environmental observations near their schools and report the data through the Internet. The GLOBE program is connected to NOAA, the National Atmospheric and Oceanic Administration.

Lab 3

Testing Water Quality

Probeware Activity

One way of judging water quality is to determine the amount of oxygen dissolved in the water. Oxygen may be supplied to a body of water from the air and from photosynthetic organisms living in the water. Clean water usually has a high oxygen content. Polluted water usually has a low oxygen content because organisms in the water use the oxygen as they decompose.

Objectives

- Using a dissolved oxygen probe, measure the concentration of dissolved oxygen in water samples obtained from different locations.
- Give reasons why the water samples have different concentrations of dissolved oxygen.

Materials

- | | | |
|--|---|--|
| <input type="checkbox"/> LabPro or CBL 2 unit | <input type="checkbox"/> Beral pipette | <input type="checkbox"/> plastic or rubber gloves |
| <input type="checkbox"/> AC adapter (optional) | <input type="checkbox"/> dissolved oxygen calibration bottle | <input type="checkbox"/> 100-mL graduated cylinder |
| <input type="checkbox"/> TI graphing calculator | <input type="checkbox"/> metric ruler | <input type="checkbox"/> 10-mL graduated cylinder |
| <input type="checkbox"/> link cable | <input type="checkbox"/> classroom thermometer | <input type="checkbox"/> 100-mL beakers (5) |
| <input type="checkbox"/> Vernier dissolved oxygen probe | <input type="checkbox"/> classroom barometer | <input type="checkbox"/> distilled water |
| <input type="checkbox"/> sodium sulfite calibration solution | <input type="checkbox"/> water samples from different locations (4 or more) | <input type="checkbox"/> wax marking pencil |
| <input type="checkbox"/> D.O. electrode filling solution | <input type="checkbox"/> jars with lids (4 or more) | <input type="checkbox"/> lab wipes |
| | | <input type="checkbox"/> laboratory apron |
| | | <input type="checkbox"/> goggles |

Procedure

Part A. Set up the Dissolved Oxygen Probe

1. Connect the TI graphing calculator to the LabPro or CBL 2 interface using the link cable. Connect the dissolved oxygen probe into Channel 1 of the interface. If the dissolved oxygen probe needs to be warmed up, proceed to Step 2. If the probe has already been warmed up, proceed to Part B.
2. Unscrew the membrane cap (counterclockwise) from the tip of the electrode on the dissolved oxygen probe. Do not touch the membrane at the very tip of the probe.
3. Use a Beral pipet to fill the membrane cap with about 1 mL of D.O. electrode filling solution. Carefully thread the membrane cap (clockwise) onto the electrode body. Do not over tighten the cap. Rinse the electrode with distilled water and carefully wipe it dry with a lab wipe.
4. Place the dissolved oxygen probe in a 250-mL beaker containing about 75 mL of water.
5. Turn on the calculator and start the DATA-MATE program. Press **CLEAR** to reset the program.
 - a. If the calculator screen displays CH 1 DO (MG/L), proceed to Step 6. If it does not, continue with this step to manually select the dissolved oxygen probe.
 - b. Select SETUP from the main screen.
 - c. Press **ENTER** to select CH 1.
 - d. Select D. OXYGEN (MG/L) from the SELECT SENSOR menu.
 - e. Select OK to return to the main screen.

Lab

3

Testing Water Quality, *continued*

Probeware Activity

6. Warm up the dissolved oxygen probe for 10 minutes.
 - a. With the probe still in the water, wait 10 minutes while the probe warms up. The probe must stay connected to the interface at all times to keep it warmed up. If disconnected for a period longer than a few minutes, it will be necessary to warm it up again.
 - b. At the end of class, leave the dissolved oxygen probe connected to the interface, with the DATAMATE program running. If this is done, the probe will stay warm and ready for the next class.




Figure 1

Part B. Calibrate the Dissolved Oxygen Probe

1. Select SETUP from the main screen.
 2. Select CALIBRATE from the setup screen.
 3. Select CALIBRATE NOW.
 4. Determine the zero-oxygen calibration point.
 - a. Remove the probe from the water and place the tip of the probe into the sodium sulfite calibration solution. **IMPORTANT:** *No air bubbles can be trapped below the tip of the probe or the probe will sense an inaccurate dissolved oxygen level.* If the voltage does not rapidly decrease, tap the side of the bottle with the probe to dislodge any bubbles. The readings should be in the 0.2- to 0.5-V range.
 - b. When the voltage stabilizes (~1 minute), press .
 - c. Enter “0” as the known concentration value in mg/L.
5. Determine the saturated DO calibration point.
 - a. Rinse the probe with distilled water and gently blot dry.
 - b. Unscrew the lid of the calibration bottle provided with the probe. Slide the lid and the grommet about 2 cm onto the probe body.
 - c. Add water to the bottle to a depth of about 1 cm and screw the bottle into the cap, as shown. **IMPORTANT:** *Do not touch the membrane or get it wet during this step.*
 - d. Keep the probe in this position for about a minute. The readings should be above 2.0 V. When the voltage stabilizes, press .
 - e. Enter the correct saturated dissolved-oxygen value (in mg/L), from the Appendix on page 41, (for example, “8.66”) using the current barometric pressure and air temperature values.
 - f. Select OK to return to the setup screen.
 - g. Select OK again to return to the main screen.
 - h. Return the dissolved oxygen probe to the beaker of water.

Part C. Finding the Dissolved Oxygen Concentration of Various Water Samples

NOTE: *Proceed with this part of the lab only after the dissolved oxygen probe has been warmed up and calibrated.*

1. In jars, collect four or more water samples from different locations. Samples could come from a tap, a pond, a lake, a river, a puddle, or an aquarium. Try to find water that has been standing and has some algae growth. Fill the jars to the top, label by source, and seal with lids. **CAUTION:** *Wear protective gloves while collecting and handling water samples.* Record your observations of the water samples in Table 1. Indicate whether any look polluted or dirty.

Lab

3

Testing Water Quality, *continued*

Probeware Activity

- With the water samples at room temperature, gently pour 25 mL of each into separate 100-mL beakers labeled with each source. Pour slowly to avoid making bubbles.
- Set up the calculator for data collection. Select SETUP from the main screen. Select MODE by pressing $\square \blacktriangle \square$ once and then pressing $\square \text{ENTER} \square$. Select SINGLE POINT from the SELECT MODE menu. Select OK to return to the main screen.
- Using a gentle motion, stir the dissolved oxygen probe through the water in one of the beakers. Make sure no bubbles are trapped under the tip of the probe. To provide an accurate reading, liquid must be continually moving past the membrane of the electrode. Once the reading displayed on the calculator screen has stabilized, select START to collect data. When data collection finishes, the dissolved oxygen concentration of the sample will be displayed on the screen.
- Record the concentration in Table 1. Press $\square \text{ENTER} \square$ to return to the main screen. Rinse the end of the probe with distilled water and place it in the next beaker to be tested.
- Repeat Step 4 for the other water samples.
- When finished, place the probe in a beaker of distilled water. Leave the DATAMATE program and calculator running for the next class. If you are the last class to use the equipment, exit the DATAMATE program and turn off the calculator. Disconnect the probe from the LabPro or CBL 2. Remove the membrane cap and rinse the inside and outside of the cap with distilled water. Rinse and carefully dry the exposed cathode and anode inner elements of the probe. Reinstall the membrane cap loosely onto the electrode body for storage.
- At the conclusion of the lab, wash your hands thoroughly with soap and water.

Data and Observations
Table 1

Sample	Water Source	Observations of Water	Concentration of Dissolved Oxygen (ppm)
1			
2			
3			
4			

Lab

3

Testing Water Quality, *continued*

Probeware Activity

Analysis

1. Explain why the water samples you collected have different concentrations of dissolved oxygen.

2. A lake sample having less than 4 ppm of dissolved oxygen is harmful to aquatic life.


a. Which of your samples could not support aquatic life?

b. Explain why oxygen dissolved in water is important for aquatic life.

3. The graph in Figure 2 shows the values for dissolved oxygen in a lake at various depths. Explain what might cause the differences in the concentrations of dissolved oxygen.

Figure 2

Oxygen Concentration at Various Depths



4. List errors you may have made in Part C that could have affected your results.

Further Investigations

- Put an aquatic plant in one water sample for a day and place the sample near a window. Determine the concentration of dissolved oxygen in the water at the beginning and end of the day, using the procedure in this lab. Explain any changes in the concentration of dissolved oxygen.
- Place the water sample containing the plant in the dark for 48 hours. Find the concentration of dissolved oxygen at the beginning and end of this period. Explain any changes.

Lab 4

How well does yeast ferment different sugars?

Probeware Activity

The most commonly used yeast, *Saccharomyces cerevisiae*, is used as baker's yeast. The function of yeast in baking is to ferment sugars found in the flour or added to the dough. This fermentation gives off carbon dioxide and ethanol. The carbon dioxide is trapped in the dough, forming gas bubbles. As the gas in the bubbles expands, the dough stretches, or rises. The baking process kills the yeast and causes the ethanol to evaporate from the dough.

In this lab, you will design an experiment to test how well yeast ferments different sugars. The amount of carbon dioxide given off by the yeast is proportional to the amount of sugar that is fermented by the yeast. The amount of carbon dioxide produced can be measured with a carbon dioxide gas sensor.

Problem

Does yeast ferment different sugars equally well, or does yeast ferment some sugars better than others?

Hypothesis

Write a hypothesis about whether yeast ferments some sugars better than others.

Objectives

- Hypothesize whether yeast ferments some sugars better than others.
- Using a carbon dioxide gas sensor, measure the amount of carbon dioxide gas produced by the fermentation of different sugars by yeast.
- Compare how well yeast ferments different sugars.

Possible Materials

- | | | |
|---|---|---|
| <input type="checkbox"/> LabPro or CBL 2 unit | <input type="checkbox"/> yeast suspension | <input type="checkbox"/> container for fermentation chamber |
| <input type="checkbox"/> carbon dioxide sensor | <input type="checkbox"/> test tube | <input type="checkbox"/> clock or watch with second hand |
| <input type="checkbox"/> TI graphing calculator | <input type="checkbox"/> stirring rod | <input type="checkbox"/> thermal mitts |
| <input type="checkbox"/> link cable | <input type="checkbox"/> 10-mL graduated cylinder | <input type="checkbox"/> laboratory apron |
| <input type="checkbox"/> AC adapter (optional) | <input type="checkbox"/> thermometer | <input type="checkbox"/> goggles |
| <input type="checkbox"/> different kinds of sugar solutions, such as glucose, sucrose (table sugar), fructose, dextrose, maltose, lactose | <input type="checkbox"/> water bath | |
| | <input type="checkbox"/> hot and cold water | |

Lab

4

How well does yeast ferment different sugars?, *continued*

Probeware Activity

Plan the Experiment

1. Decide on a procedure that you can use to compare how well yeast ferments different sugars. The procedure should include a probeware system consisting of a carbon dioxide gas sensor, LabPro or CBL 2 interface, graphing calculator, and link cable. Other possible materials are listed on the previous page.
2. Prepare a water bath to incubate the yeast suspension. The temperature of the water bath should be maintained at a constant temperature of around 37°C.
3. Design a fermentation chamber in which the yeast will ferment the sugar.
4. Decide which sugars you will test. Think about how much yeast suspension and sugar solution you will use. After mixing the yeast suspension with the sugar solution, allow the mixture to incubate at 37°C for 10 minutes before adding it to the fermentation chamber.
5. Observe the mixture in the fermentation chamber. Use caution when placing the carbon dioxide probe into the fermentation chamber. Twist the stopper slightly to seal it in the opening being careful not to twist the shaft of the sensor. Use the carbon dioxide sensor to collect data on the amount of carbon dioxide gas produced from the fermentation of the sugar by the yeast.
6. To collect data, start the DATAMATE program and go to the MAIN MENU. Press **CLEAR** to reset the program. If the DATAMATE program is not loaded, transfer the program from the memory of the LabPro or CBL 2 interface to the TI graphing calculator. Do not select START until you are ready to collect data from the sensor in the fermentation chamber.

7. Decide how frequently you will collect data and for how long. You can record your observations and data in the tables provided or make your own tables. Label the columns appropriately.
8. Write your procedure on another sheet of paper or in your notebook. It should include the amounts of each material you will need.




Figure 1

Check the Plan

1. Be sure that you include a control fermentation chamber in your experiment and that the experimental fermentation chambers vary in one way only.
2. Make sure your teacher has approved your experimental plan before you proceed further.
3. Carry out the experiment.
4. As you do the experiment, keep the water bath at 37°C.
5. Completely rinse the fermentation chamber with water after each trial to remove used sugar, yeast, and carbon dioxide. Thoroughly dry the inside with a paper towel. Fan the carbon dioxide sensor with air after each trial. This will clear the sensor of any excess carbon dioxide.

Lab

4

How well does yeast ferment different sugars?, *continued***Probeware Activity****Data and Observations****Table 1**

Type of Sugar in Chamber	Observations

Table 2

Carbon Dioxide Concentration (ppm)				
Time (in seconds)	Sugar used: _____	Sugar used: _____	Sugar used: _____	Sugar used: _____

Table 3

Type of Sugar Used	Rate of Fermentation (ppm/s)

Lab

4

How well does yeast ferment different sugars?, *continued***Probeware Activity****Analysis**

1. Use the graphing calculator to graph carbon dioxide concentration (ppm) versus time for the control trial and for each sugar trial.
 - a. Select ANALYZE from the main screen.
 - b. Select CURVE FIT from the ANALYZE OPTIONS menu.
 - c. Select LINEAR (CH 1 VS TIME) from the CURVE FIT menu.
 - d. The linear-regression statistics for these two lists are displayed for the equation in the form:
$$Y = A \cdot X + B$$
 - e. Enter the value of the slope, A , as the rate of fermentation in Table 3.
 - f. Press **ENTER** to view a graph of the data and the regression line.
 - g. Press **ENTER** to return to the ANALYZE menu.
 - h. Select RETURN TO MAIN SCREEN from the ANALYZE menu. Repeat the procedure for each trial.

2. How well did the yeast ferment the different sugars?

3. What is a possible explanation for your answer to question 2?

4. Why was it important to keep the water bath in which the yeast was incubated at 37°C?

5. What variables did you keep constant in your experiment?

6. Explain the data you recorded for the control fermentation chamber.

Checking Your Hypothesis

Was your hypothesis supported by your data? Why or why not?

Further Investigations


1. Repeat the lab using different types of baker's yeast, such as dry yeast, cake yeast, and quick-starting yeast. Compare the fermentation rates of the different yeasts.
2. Research the genome map of yeast and present this information to the class.

Lab 5

How can pH be used to compare rates of photosynthesis?

Probeware Activity

During respiration, aquatic plants, such as *Elodea*, release carbon dioxide into the water. During photosynthesis, the plants extract the carbon dioxide they need from the water. Carbon dioxide dissolves in water to form a weak acid called carbonic acid (H_2CO_3). Some of the carbonic acid dissociates to form H^+ ions and HCO_3^- ions, as shown below.



The double direction arrows indicate that the chemical reactions can proceed in either direction, and thus are reversible. As the carbon dioxide concentration in water decreases, the concentration of carbonic acid also decreases, which increases the pH of the water. Conversely, an increase in the water's carbon dioxide concentration results in a decrease in pH. Therefore, the pH of water can be used as an indication of the level of photosynthesis of aquatic plants. The greater the rate of photosynthesis, the greater the amount of dissolved carbon dioxide used, and thus the higher the pH of the water. In this lab, you will design an experiment to test the effect of light intensity on the photosynthetic rate of an aquatic plant, using a probe that measures pH.

Problem

Do plants undergo different rates of photosynthesis at different times of the day?
Do plants undergo different rates of photosynthesis in different parts of the ocean? How does light intensity affect the rate of photosynthesis of plants?

Hypothesis

Write a hypothesis about the effect of light intensity on the rate of photosynthesis of aquatic plants.

Objectives

- Make a hypothesis about the effect of light intensity on the rate of photosynthesis.
- Expose an aquatic plant to different intensities of light.
- Use a probe to measure the pH of water samples.
- Compare the rates of photosynthesis of an aquatic plant exposed to different light intensities.

Lab

5

How can pH be used to compare rates of photosynthesis?, *continued*

Probeware Activity

Possible Materials



- | | | |
|---|---|--|
| <input type="checkbox"/> LabPro or CBL 2 unit | <input type="checkbox"/> aquatic plant, such as <i>Elodea</i> | <input type="checkbox"/> rinse bottle of distilled water |
| <input type="checkbox"/> TI graphing calculator | <input type="checkbox"/> large test tubes with stoppers | <input type="checkbox"/> 250-mL beaker |
| <input type="checkbox"/> link cable | <input type="checkbox"/> test-tube rack | <input type="checkbox"/> laboratory apron |
| <input type="checkbox"/> AC adapter (optional) | <input type="checkbox"/> various light sources | <input type="checkbox"/> goggles |
| <input type="checkbox"/> pH sensor | <input type="checkbox"/> dechlorinated water | |

Plan The Experiment

- Decide on a procedure that you can use to compare the rates of photosynthesis of an aquatic plant exposed to different intensities of light. The procedure should use a probeware system that includes a pH probe, LabPro or CBL 2 interface, TI graphing calculator, and link cable. You may wish to use the other suggested materials as well.
- Think about how you will vary light intensity. How many samples will you test? Decide on the variables that you will need to keep constant during the experiment. What will be your control? **CAUTION: To avoid harming the aquatic plants you will be working with, do not place the samples too close to an artificial light source. Heat from the light source may increase the temperature of the water above the plant's tolerance level.**
- Decide how frequently you will collect data and for how long. To collect data, plug the pH sensor into Channel 1 of the LabPro or CBL 2 interface. Using the link cable, connect the TI graphing calculator to the interface. Push the link cable securely into each jack.
- Turn on the graphing calculator. Start the DATAMATE program and go to the MAIN MENU. Press **CLEAR** to reset the program. If the DATAMATE program is not loaded, transfer the program from the memory of the LabPro or CBL 2 interface to the TI graphing calculator.
- Set up the calculator and interface for a pH sensor.
 - If the calculator displays PH in CH 1, proceed directly to Step 6. If it does not, continue with this step to set up your sensor manually.
 - Select SETUP from the main screen.
 - Press **ENTER** to select CH 1.
 - Select PH from the SELECT SENSOR menu.
 - Select OK to return to the main screen. Readings from the pH sensor will be displayed on the main screen.
- Before each use of the pH probe, rinse the tip of the electrode completely with distilled water. Carefully hold the pH probe over a beaker and use a rinse bottle to gently rinse the tip with distilled water. **CAUTION: The tip of the pH probe is fragile and can be broken easily.** Do not let the probe dry out. When not in use during the lab, keep the probe immersed in a beaker of tap water.
- Place the probe in the sample. Allow the pH value reading to stabilize for 15 to 20 seconds. Record the value in Table 1, or you can make your own data table.
- Write your procedure on another sheet of paper or in your notebook. It should include all the materials you will use.

Lab

5

How can pH be used to compare rates of photosynthesis?, *continued*

Probeware Activity




Figure 1

Check The Plan

1. Be sure that light intensity is the only variable that changes in your experiment.
2. Did you measure the pH of your samples before beginning the experiment?
3. Make sure the teacher has approved your experimental plan before you proceed further.
4. Carry out your experiment.
5. After the experiment, return the plants to their original location or as directed by the teacher.

Data and Observations

Table 1

pH of Samples			
Day	Type of light intensity: _____	Type of light intensity: _____	Type of light intensity: _____

Analysis

1. What variables did you keep constant in your experiment?

2. Which sample had the highest pH? The lowest pH?

Lab

5**How can pH be used to compare rates of photosynthesis?, *continued*****Probeware Activity**

3. Explain the differences in the pH values of the samples.

4. During what time of day would you expect outdoor plants to undergo the highest rate of photosynthesis? Explain.

5. In what part of the ocean would you expect to find the most photosynthetic organisms? Explain.

Checking Your Hypothesis

Was your hypothesis supported by your data? Why or why not?

Further Investigations

1. Repeat the experiment to investigate how the color of light affects the rate of photosynthesis.
2. Repeat the experiment devising a method to investigate whether temperature has an effect on the rate of photosynthesis.

Lab 6

Effect of Environmental Temperature on the Metabolic Rates of Animals

Probeware Activity

Endothermic (warm-blooded) animals, such as mammals, have the ability to maintain a constant internal body temperature through metabolic activity, but ectothermic (cold-blooded) animals do not. They do not generate much heat through their metabolism, so they must pick up heat from their environments. An ectothermic animal controls its internal temperature through its behavior. It may bask in the sun to get warmer or burrow underground to cool off. Ectothermic animals are more active in warm temperatures, and they slow down when the environment becomes cold. The rate of their activity can be examined by determining their metabolic rate. Animals consume food and break it down, producing carbon dioxide in the process of respiration. The overall reaction can be summarized by the following equation:



The rate of carbon dioxide production can be used as a measure of the metabolic rate.

In this lab, you will be studying the relationship between the temperature of the environment and the metabolic rate of ectothermic animals. You will be able to monitor metabolic rate by measuring the carbon dioxide production during respiration, using a carbon dioxide gas sensor.

Objectives

- Use a carbon dioxide gas sensor to measure carbon dioxide concentrations.
- Determine the metabolic rate of small animals by measuring their rate of carbon dioxide production.
- Explain the effect of temperature on the metabolic rate of cold-blooded animals.

Materials

- | | | |
|---|---|--|
| <input type="checkbox"/> LabPro or CBL 2 unit | <input type="checkbox"/> thermometer | <input type="checkbox"/> Beral pipette or small beaker |
| <input type="checkbox"/> TI graphing calculator | <input type="checkbox"/> crickets or pill bugs (4–10) | <input type="checkbox"/> balance |
| <input type="checkbox"/> link cable | <input type="checkbox"/> warm tap water | <input type="checkbox"/> paper towels |
| <input type="checkbox"/> carbon dioxide sensor | <input type="checkbox"/> ice | <input type="checkbox"/> laboratory apron |
| <input type="checkbox"/> 250-mL respiration chamber (comes with sensor) | <input type="checkbox"/> 500-mL or 1000-mL beaker | <input type="checkbox"/> goggles |

Lab

6

Effect of Environmental Temperature on the Metabolic Rates of Animals, *continued*

Probeware Activity

Procedure

To divide up the work for this lab, lab teams can work in three groups to test the same type of animal at different temperatures. Group One should collect data using temperatures at intervals between 10 and 20°C. Group Two should collect data using temperatures at intervals between 20 and 30°C. Group Three should collect data using temperatures at intervals between 30 and 40°C. Although the three groups will be using different individual animals that may differ in mass, the slight differences in mass of the animals should not make a significant difference in the results. An individual animal with greater mass should produce carbon dioxide at a higher rate, but we will assume that the rate will be proportional to its difference in mass. The carbon dioxide production rate, measured in ppm of CO₂/gram of animal mass, will be computed. At the conclusion of the lab, the three lab groups should share data.

1. Connect the LabPro or CBL 2 unit and TI graphing calculator with a link cable. Be sure to press the cable firmly into each piece of equipment. Connect the carbon dioxide sensor into CH 1 of the LabPro or CBL 2.
2. Turn on the graphing calculator and start the DATAMATE program. Press **CLEAR** to reset the program.
3. Set up the calculator and interface for a CO₂ gas sensor.
 - a. Select SETUP from the main screen.
 - b. If the calculator displays CO2 GAS (PPM) in CH 1, proceed directly to Step 4. If it does not, continue with this step to set up your sensor manually.
 - c. Press **ENTER** to select CH 1.
 - d. Select CO2 GAS from the SELECT SENSOR menu.
 - e. Select parts per million (PPM) as the unit.
4. Set up the data collection parameters.
 - a. Press **▲** once and then press **ENTER** to select MODE.
 - b. Select TIME GRAPH from the MODE menu.
 - c. Select CHANGE TIME SETTINGS from the TIME GRAPH SETTINGS menu.
 - d. Enter “10” as the time between samples in seconds.
 - e. Enter “60” as the number of samples to be collected.
 - f. Select OK to return to the setup screen.
 - g. Select OK to return to the main screen.
5. Prepare a water bath according to your assigned temperature.
 - **Group 1:** To prepare a 10 to 20°C water bath, place cool tap water in the 1000-mL beaker, add ice, and mix. Measure the temperature of the mixture. Add ice or cool water and mix until the water reaches the desired temperature.
 - **Group 2:** To prepare a 20 to 30°C water bath, start with warmer tap water and cool it with ice as needed.
 - **Group 3:** To prepare a 30 to 40°C water bath, start with warmer tap water and warm it with hot water as needed.
6. Measure the mass of your respiration chamber. Place your animals in the respiration chamber, and measure the mass of the respiration chamber again. **CAUTION: Handle live animals with care.** Calculate the difference between your measurements to find the mass of your animals. Record the mass of the animals in Table 1.
7. Carefully put the carbon dioxide probe into the top of the respiration chamber. Twist the stopper slightly to seat it in the opening. Be careful not to twist the shaft of the sensor. Place the respiration chamber in the water bath as shown in Figure 1. This setup will keep the animals at a constant temperature throughout the lab.
8. Wait five minutes to allow the temperature of the respiration chamber to adjust to the water bath temperature. Measure and record the water bath temperature in Table 1.

Lab

6

Effect of Environmental Temperature on the Metabolic Rates of Animals, *continued*

Probeware Activity

Select START from the main screen to begin data collection. While data is being collected, try to maintain a constant temperature in your water bath.

9. When data collection has finished, a graph of CO₂ GAS vs. TIME will be displayed. Press **ENTER** to return to the main screen.
10. Calculate the rate of respiration.
 - a. Select ANALYZE from the main screen.
 - b. Select CURVE FIT from the ANALYZE OPTIONS menu.
 - c. Select LINEAR (CH 1 VS TIME) from the CURVE FIT menu.
 - d. The best-fit line for your data is displayed in the form: $Y = A \cdot X + B$
- e. Enter the value of the slope, A , as the rate of carbon dioxide production in your data table.
- f. Press **ENTER** to view a graph of the data and the regression line.
- g. Press **ENTER** to return to the ANALYZE menu.
- h. Select RETURN TO MAIN SCREEN from the ANALYZE menu.
11. Remove the respiration chamber from the water bath. Dry off the respiration chamber and place it on the lab table. Carefully remove the carbon dioxide probe from the respiration chamber.
12. After all the data has been collected, place your animals in the location designated by your teacher. Exit the DATAMATE program and turn off the calculator. Share your data with the other lab groups.




Figure 1

Data and Observations

Animal Type _____

Table 1

Group	Animal Mass (g)	Temperature (°C)	CO ₂ Production Rate (ppm/s)	CO ₂ Production Rate per gram (ppm/s · g)
1				
2				
3				

Lab

6

Effect of Environmental Temperature on the Metabolic Rates of Animals, *continued***Probeware Activity****Analysis**

1. Determine the rate of CO₂ produced per gram of animal mass. Divide the CO₂ production rate by the animal mass. Record the CO₂ production rate per gram in Table 1.

2. Plot a graph of CO₂ production rate per gram versus temperature, plotting the CO₂ production rate (ppm/s · g) on the *y*-axis and temperature (°C) on the *x*-axis.

3. Compare the carbon dioxide production rates at the various temperatures tested. At what temperatures was the production rate of carbon dioxide the lowest and the highest?

4. What can you infer about the effect of temperature on metabolic rate?

5. Lizards are ectothermic animals. Suggest a possible explanation for a lizard being able to change its skin color other than for camouflage as protection against predators.

Further Explorations

1. Repeat the experiment with different types of animals such as cockroaches, flies, pill bugs, or mealworms. Compare the respiration rates for these animals.
2. Do research to find out more about the ways that ectothermic animals control their internal body temperature through behavior.

**Lab
7****What is the effect of exercise on body temperature?****Probeware Activity**

The normal body temperature of a human is about 37°C. Normal metabolic processes generate heat that warms the body. The rate of heat production can be increased or decreased by a change in muscle activity or by hormonal action that changes metabolic rate. Sweating and the dilation of blood vessels help to cool the body when it becomes overheated. Shivering and the constriction of blood vessels help to warm the body when it becomes too cold. In this lab, you will design an experiment to test the effect of exercise on a person's body surface temperature. To monitor the temperature, you will attach a temperature probe to the skin.

Problem

How does exercise affect body surface temperature?

Hypothesis

Hypothesize how exercise affects a person's body surface temperature. Write your hypothesis on the lines below.

Objectives

- Hypothesize about the effect of exercise on body surface temperature.
- Measure body surface temperature using a temperature probe.
- Compare body surface temperature while at rest with body surface temperature while exercising.

Materials   

- | | | |
|---|--|---|
| <input type="checkbox"/> LabPro or CBL 2 unit | <input type="checkbox"/> athletic tape | <input type="checkbox"/> plastic bag to cover the temperature probe |
| <input type="checkbox"/> TI graphing calculator | <input type="checkbox"/> material for insulating the temperature probe, such as wool or polyester fill | <input type="checkbox"/> laboratory apron |
| <input type="checkbox"/> link cable | | <input type="checkbox"/> goggles |
| <input type="checkbox"/> temperature probe | | |

Lab

7

What is the effect of exercise on body temperature?, *continued*

Probeware Activity

Plan the Experiment

1. Decide on a procedure that you can use to test the effects of exercise on body surface temperature. The procedure should use a probeware system that includes a temperature probe, LabPro or CBL 2 interface, and TI graphing calculator.
2. As you develop your procedure, think about where on the body you will measure surface temperature. You might measure the surface temperature in different locations on the body before choosing a location to monitor while exercising. You can use Table 1 to record your resting-temperature data, or make your own table.
3. Choose an exercise that is safe for you to do.
CAUTION: *If you have any conditions that may be aggravated by exercise, inform the teacher.*
4. Cover the temperature probe with a plastic bag to prevent sweat from getting on the probe.
5. Think about how often you will collect data while exercising and for how long. What will be your control? To collect data, set up the probeware system by connecting the link cable between the graphing calculator and the LabPro or CBL 2 interface. Firmly press the cable into each unit. Plug the temperature probe into Channel 1 of the interface. Turn on the calculator and start the DATAMATE program.
6. Place the probe against the skin in the selected location. If part of the probe is exposed to the air, cover that part with an insulating material.
7. Temperature can be continuously monitored or collected for a specific amount of time. The current temperature readings are displayed on the calculator screen and updated once a second. If you wish to collect data for a period of time, select START. The DATAMATE program will collect temperature readings for 180 seconds. If you wish to collect data for a longer time, select SETUP

from the main screen. Press \uparrow to select MODE and press ENTER . Select TIME GRAPH from the MODE menu. Select CHANGE TIME SETTINGS from the TIME GRAPH menu. Enter how often you want data collected and how many points you wish to collect. To return to the main screen, select OK from the TIME GRAPH and SETUP menu. Select START when you want to begin data collection.

8. When data collection has finished, the calculator will display the graph of temperature versus time. After you are done looking at the graph, press ENTER to return to the main screen. Exit the program and press STAT then select EDIT to view your data. The time values will be in list L1 and the temperature readings will be in list L2. You can record selected data in Table 2, or make your own table.
9. Write your procedure on another sheet of paper or in your notebook. It should include any types of exercise equipment you will use.




Figure 1

Check the Plan

1. Your experiment should include a control.
2. Choose an exercise that is safe. Be sure that the exercise causes little disruption to other groups in the classroom.
3. Make sure the teacher has approved your experimental plan before you proceed further.
4. Carry out your experiment.

Lab

7

What is the effect of exercise on body temperature?, *continued***Probeware Activity****Analysis**

1. Describe the control in your experiment. What was its purpose?

2. Why might some parts of the body have a higher surface temperature than others?

3. What is the effect of exercise on body surface temperature? Why do you think exercise has this effect?

4. Explain the effect of the body's cooling mechanisms, such as blood-vessel dilation and sweating, on the data you obtained.

5. Would your data have differed if you had worn insulating clothing while exercising? Explain.

Checking Your Hypothesis

Was your hypothesis supported by your data? Why or why not?

Further Investigations

1. Repeat the experiment. This time collect body surface temperature data over a longer period of time and monitor the skin for sweating and a flushed appearance. Compare the temperature data you collect with those obtained in the first experiment.
2. Test the effects of different exercises on body surface temperature. Do some exercises affect body temperature more than others? If so, why?

Lab
8

Measuring Response Time

Probeware Activity

If you were walking down the street and heard someone call your name, it would take you a short time to look up, determine where the sound came from, and respond by waving. This interval of time is called your response time. You need time to perceive and process any stimulus, whether it is something you see, smell, feel, taste, or hear. For this lab, you will determine the time it takes you to respond to a sound stimulus.

Objectives

- Use a motion detector to determine the time it takes for a person to respond after hearing a sound.
- Conduct repeated trials and calculate an average response time.
- Compare the response times of different individuals.

Materials   

- C-clamp
- LabPro or CBL 2 unit
- TI graphing calculator
- link cable
- motion detector

Procedure

1. Connect the LabPro or CBL 2 unit to the TI graphing calculator, using the link cable. Connect the motion detector to the DIG/SONIC 1 port on the right side of the interface. Place the motion detector about 4 meters from a blank wall, facing the wall, as shown in Figure 1.
2. Turn on the calculator and start the DATAMATE program. Press **CLEAR** to reset the program. If the program has not been loaded on the graphing calculator, transfer the program from the memory of the LabPro or CBL 2.
3. One of the students in the group will serve as the runner. That student should stand at the position shown in Figure 1, about 2 meters from the motion detector and slightly to one side. When data collection begins, the motion detector will transmit ultrasonic sound waves toward the wall. The runner should start outside of the range of the beam.

Lab

8

Measuring Response Time, *continued*

Probeware Activity




Figure 1




Figure 2

- Another student will serve as the operator of the calculator and motion detector. The operator will select START to begin collecting data for the activity.
- The motion detector makes a clicking sound once it begins emitting sound waves. As soon as the runner hears the clicking sound, he or she will run across the beam of the motion detector, as shown in Figure 1. Since the experiment seeks to measure the runner's response time to the sound, the runner should not run in response to visual stimuli. For that reason, the runner should not look at the operator or at the motion detector, since some models have an LED that becomes lit when the motion detector is on.
- For 5 seconds after the motion detector has been turned on, the calculator will collect data. When data collection has finished, press **ENTER** to display a graph of distance vs. time. The graph may look like that shown in Figure 2. The horizontal part of the line indicates that the motion detector is bouncing sound waves off the wall, a stationary object. The depression in the line indicates that the sound waves have bounced off something closer to the detector—the runner as he or she ran across the beam.
- Use the arrow keys on the calculator to trace along the curve until it gets to the lowest part of the depression. This point represents the time at which the runner crossed the beam. The x -value and y -value for that point will show on the screen. The x -value is the response time, or the time it took for the student to respond to the clicking of the motion detector and move across the beam. The y -value at this point indicates the distance between the motion detector and the runner. Record the x -value in Table 1 of Data and Observations, in the space for Trial 1.
- When you are done working with the graph on the graphing calculator, press **ENTER**. Select MAIN SCREEN from the graph menu and repeat Steps 3–8 to collect more data. Every student in the group should be tested, and each student should do at least 5 trials. Make sure each student starts in the same position and that the motion detector is held in the same position for all of the trials. Copy the graph for one of your trials onto the grid in Data and Observations. (If you can print the LCD screen using TI-GRAPH LINK™, you can do that instead.)
- When you have finished collecting data for all students in the group, select QUIT from the main screen and turn off the graphing calculator.

Analysis

1. Which student had the fastest response time in your group?

2. Did response times improve as more trials were done? Explain why you think improvement did or did not occur.

3. Was there any difference in the response time of students of different ages?

4. How do you think distance from the source of a sound would affect response time?

5. The speed of sound in air is 343 meters per second at 20°C. The speed of sound increases as air temperature increases at a rate of 0.6 meters per second for each degree Celsius increase. Do you think an increase in air temperature would significantly increase your response time in this experiment? Why or why not?

Further Explorations

1. Use this experiment to compare response times between students and teachers.
2. Do the same experiment, but use the LED light on the ultrasonic motion detector as a visual stimulus to start running. Compare visual response time to auditory response time.

Lab 9

Breathing and Heart Rate

Probeware Activity

Was there ever a time when you were under a condition of high emotional stress and started to breathe rapidly? An excessive increase in breathing rate is called **hyperventilation**. A decrease in breathing rate is called **hypoventilation**. In this lab, you will examine how changes in breathing rate can affect a person's heart rate, or the number of times the heart beats per minute.

Objectives

- Use a heart rate monitor to measure your heart rate.
- Determine how changes in your breathing rate affect your heart rate.

Materials

- | | | |
|--|--|---|
| <input type="checkbox"/> LabPro or CBL 2 unit | <input type="checkbox"/> link cable | <input type="checkbox"/> saline solution |
| <input type="checkbox"/> Vernier Exercise Heart Rate Monitor | <input type="checkbox"/> stopwatch or clock with second hand | <input type="checkbox"/> colored pencils (2 colors) |
| <input type="checkbox"/> TI graphing calculator | <input type="checkbox"/> chair | <input type="checkbox"/> laboratory apron |
| | | <input type="checkbox"/> goggles |

Procedure

Part A. Hyperventilation

1. Connect the LabPro or CBL 2 unit to the graphing calculator using a link cable. Connect the receiver module of the exercise heart rate monitor to Channel 1. Turn on the graphing calculator and start the DATAMATE program. Press **[CLEAR]** to reset the program.
2. Elastic straps, for securing the transmitter belt, come in two different sizes (small and large). Select the size of elastic strap that best fits the subject being tested. It is important that the strap provides a snug fit of the transmitter belt. Wet each of the electrodes (the two grooved rectangular areas on the underside of the transmitter belt) with 3 drops of saline solution. Secure the transmitter belt against the skin directly over the base of the rib cage. The POLAR logo on the front of the belt should be centered on the chest. Adjust the elastic strap to ensure a tight fit.
3. Set up the calculator and interface for an exercise heart rate monitor.

- a. If the calculator displays HEART RT (BPM) in CH 1, proceed directly to Step 4. If it does not, continue with this step to set up your sensor manually.
- b. Select SETUP from the main screen.
- c. Press **[ENTER]** select CH 1.
- d. Select HEART RATE from the SELECT SENSOR menu.
- e. Select EX HEART RATE from the HEART RATE menu.
- f. Select OK to return to the main screen.




Figure 1

Lab

9

Breathing and Heart Rate, *continued*

Probeware Activity

4. Determine that the exercise heart rate monitor is working properly. Heart rate readings displayed on the calculator screen should be steady (± 8 beats per minute) and within the normal range of the individual—usually between 60 and 100 beats per minute. If readings appear sporadic, tighten the elastic strap and move the receiver module closer to the center of the transmitter belt. If all equipment is functioning properly, proceed to Step 5.
5. Sit on a chair and breathe at a normal rate. Select START to begin collecting heart rate in beats per minute (BPM). The program displays the average heart rate every 5 seconds on the screen and plots a graph.
6. After collecting data at a normal breathing rate for 60 seconds, continue collecting data while breathing rapidly for the next 30 seconds. **CAUTION: If you feel light-headed or dizzy, return to normal breathing immediately.** Breathe normally for another 30 seconds, while continuing to collect data. Then press **STO▶** to stop data collection. A graph of heart rate vs. time will be displayed. Press **▶** to move the trace cursor from one data point to the next. The values of each data point are displayed at the bottom of the calculator screen. Record the data in Table 1. Plot your data from Table 1 on the graph in Data and Observations,

using a colored pencil. Plot time on the x -axis and BPM on the y -axis.

7. Press **ENTER** to return to the main screen.
8. Repeat Part A, steps 2–7, for the other members of your group.

Part B. Hypoventilation

1. Repeat Part A, steps 2–4.
2. After collecting data at a normal breathing rate for 60 seconds, continue collecting data while holding your breath for the next 30 seconds. **CAUTION: Take a breath if you need to before the 30 seconds are up.** Breathe normally for another 30 seconds, while continuing to collect data, and then press **STO▶** to stop data collection. A graph of Heart Rate vs. Time will be displayed.
3. Press **▶** to move the trace cursor from one data point to the next. The values of each data point are displayed at the bottom of the calculator screen. Record the data in Table 1. Plot your data from Table 1 on the graph in Data and Observations, using a different colored pencil than the one you used in Part A. Plot time on the x -axis and BPM on the y -axis.
4. Press **ENTER** to return to the main screen.
5. Repeat Part B, steps 1–4, for the other members of your group.

Data and Observations**Table 1**

Time (s)	Hyperventilation Heart Rate (BPM)	Hypoventilation Heart Rate (BPM)
15		
30		
45		
60		
75		
90		
105		
120		

Lab
9

Breathing and Heart Rate, *continued*

Probeware Activity

A large grid of graph paper, consisting of 20 columns and 30 rows of small squares, intended for recording data during the lab activity.

Lab

9

Breathing and Heart Rate, *continued***Probeware Activity****Analysis**

1. What effect does hyperventilation have on heart rate?

2. What effect does hypoventilation have on heart rate?

3. When athletes undergo strenuous exercise, they are told to try to breathe at a normal rate. Why is this so?

4. Extreme emotional stress, fright, or pain can cause hyperventilation. Why is a person who is undergoing hyperventilation told to “catch your breath”?

5. In what situations would it be a survival advantage for an animal to undergo hypoventilation?

Further Explorations

1. Your breathing pattern (sequence and duration of exhalation and inhalation) can change during the course of a day. Sleeping, coughing, yawning, sneezing, laughing, hiccuping, and crying cause changes in breathing patterns. Describe in your own words the breathing pattern you experience during each act. Compare your descriptions with those of your classmates. Discuss any differences in the descriptions.
2. Research the disease called Congenital Central Hypoventilation Syndrome, sometimes called Ondine Curse. Write a report that includes a discussion of the causes and symptoms of the disease. Present your findings to the class.

Lab 10

What is the effect of exercise on heart rate?

Probeware Activity

The heart pumps blood to all the cells of the body. As a person exercises, the muscle cells use increased amounts of oxygen and food, which must be replaced by the blood. The muscle cells also produce more wastes, which must be removed by the blood. In this lab, you will design an experiment to test the effects of exercise and other factors on a person's heart rate, which is the number of heartbeats per minute. You will use a heart rate monitor to measure heart rate.

Problem

Does a shift in body position cause a person's heart rate to change? How does exercise affect heart rate? After exercising, how long does it take for a person's heart to recover, that is, return to its resting rate?

Hypothesis

Write a hypothesis about how exercise and body position affect a person's heart rate.

Objectives

- Hypothesize about the effects of body position and exercise on heart rate.
- Use a heart rate monitor to measure heart rate.
- Compare the effects of different body positions and exercises on heart rate.
- After performing different exercises, determine the heart's recovery time.

Possible Materials

- | | | |
|---|--|---|
| <input type="checkbox"/> LabPro or CBL 2 unit | <input type="checkbox"/> Vernier Exercise Heart Rate Monitor | <input type="checkbox"/> saline solution |
| <input type="checkbox"/> TI graphing calculator | <input type="checkbox"/> stopwatch or clock with second hand | <input type="checkbox"/> exercise equipment |
| <input type="checkbox"/> link cable | | <input type="checkbox"/> laboratory apron |
| | | <input type="checkbox"/> goggles |

Lab

10

What is the effect of exercise on heart rate?, *continued*

Probeware Activity

Plan the Experiment

1. Devise a procedure that tests the effects of body position and various exercises on heart rate, and determines the heart's recovery time after each exercise. The procedure should include a probeware system consisting of an exercise heart rate monitor, LabPro or CBL 2 interface, TI graphing calculator, and link cable.
2. Choose body positions and exercises that are safe for you. **CAUTION: If you have any medical conditions that may be aggravated by these movements, inform your teacher.** In your procedure, include how long you will maintain each body position and perform each exercise, as well as how many measurements you will make. Think about what your control will be in the experiment.
3. Connect the LabPro or CBL 2 unit to the graphing calculator using a link cable. Connect the receiver module of the exercise heart rate monitor to CH 1. Turn on the graphing calculator and start the DATAMATE program. Press **CLEAR** to reset the program.
4. Elastic straps, for securing the transmitter belt, come in two different sizes (small and large). Select the size of elastic strap that best fits the subject being tested. It is important that the strap provides a snug fit of the transmitter belt. Wet each of the electrodes (the two grooved rectangular areas on the underside of the transmitter belt) with 3 drops of saline solution. Secure the transmitter belt against the skin directly over the base of the rib cage. The POLAR logo on the front of the belt should be centered on the chest. Adjust the elastic strap to ensure a tight fit.
5. Set up the calculator and interface for an Exercise Heart Rate Monitor.
 - a. If the calculator displays HEART RT (BPM) in CH 1, proceed directly to Step 6. If it does not, continue with this step to set up your sensor manually.
 - b. Select SETUP from the main screen.
 - c. Press **ENTER** select CH 1.
 - d. Select HEART RATE from the SELECT SENSOR menu.

- e. Select EX HEART RATE from the HEART RATE menu.
- f. Select OK to return to the main screen.




Figure 1

6. Determine that the exercise heart rate monitor is working properly. Heart rate readings displayed on the calculator screen should be steady (± 8 beats per minute) and within the normal range of the individual—usually between 60 and 100 beats per minute. If readings appear sporadic, tighten the elastic strap and move the receiver module closer to the center of the transmitter belt. If all equipment is functioning properly, proceed to Step a.
 - a. Plan each trial. Trials should include measuring the average resting heart rate for several different body positions.
 - b. For exercise trials, be sure to continue collecting data until the resting heart rate is reached. Recovery time is the amount of time it takes to reach the resting heart rate after exercise has stopped.
 - c. Select START to begin collecting heart rate in beats per minute (BPM). The program displays the average heart rate every 5 seconds on the screen and plots a graph.
 - d. Start the exercise to be tested. A graph of heart rate vs. time will be plotted while the exercise is being performed.
 - e. When you are finished exercising, wait until the resting heart rate is reached and press **STO▶** to stop data collection.

**Lab
10****What is the effect of exercise on heart rate?, *continued*****Probeware Activity**

A scaled graph of heart rate vs. time will be displayed. Press to move the trace cursor from one data point to the next. The values of each data point are displayed at the bottom of the calculator screen. Record the data in Tables 1 and 2.

- f. Press to return to the main screen.
- g. Repeat steps c–f for each heart rate trial.
7. Use Tables 1 and 2 to record your data, or you may wish to make your own data tables.
8. Write your procedure on another sheet of paper or in your notebook. It should include any types of exercise equipment you use.

Check the Plan

1. Your experiment should include a control.
2. Be sure that the exercises you do are not time-consuming and cause little disruption to other groups in the classroom.
3. Choose exercises that are safe.
4. Make sure your teacher has approved your experimental plan before you proceed further.
5. Carry out the experiment.

Data and Observations**Table 1**

Body Position	Average Resting Heart Rate (BPM)

Table 2

Effect of Exercise on Heart Rate			
Type of exercise	Duration of exercise (minutes)	Maximum heart rate during exercise (BPM)	Recovery time (seconds)

Lab

10

What is the effect of exercise on heart rate?, *continued***Probeware Activity****Analysis**

1. What was the control in your experiment? What did it show?

2. What body position resulted in the highest heart rate? Explain why.

3. What type of exercise increased your heart rate the most? The least?

4. Explain why the exercises you identified in your answer to question 3 had different effects on heart rate.

5. What is the relationship between type of exercise, heart rate, and recovery time?

6. Two people who have performed the same exercise for the same amount of time may show different recovery times. Why do you think this is so?

7. What factors, other than exercise and body position, may affect heart rate?

Checking Your Hypothesis

Was your hypothesis supported by your data? Why or why not?

Further Investigations

1. Pool class data to compare the average heart rate of males with that of females in your class.
2. Use the probeware system you used in the lab to observe the effects of different factors on heart rate. You might consider testing the effects of time of day, type or amount of food consumed, or stress.

Appendix

Dissolved Oxygen Concentrations

Use this table to calibrate the dissolved oxygen probe used in Lab 3 *Testing Water Quality*.

Dissolved Oxygen (mg/L) in Air-Saturated Distilled Water								
Air Temperature	Barometric Pressure							
	770 mm Hg	760 mm Hg	750 mm Hg	740 mm Hg	730 mm Hg	720 mm Hg	710 mm Hg	700 mm Hg
17°C	9.86	9.74	9.61	9.48	9.35	9.22	9.10	8.97
18°C	9.67	9.54	9.41	9.29	9.16	9.04	8.91	8.79
19°C	9.47	9.35	9.23	9.11	8.98	8.86	8.74	8.61
20°C	9.29	9.17	9.05	8.93	8.81	8.69	8.57	8.45
21°C	9.11	9.00	8.88	8.76	8.64	8.52	8.40	8.28
22°C	8.94	8.83	8.71	8.59	8.48	8.36	8.25	8.13
23°C	8.78	8.66	8.55	8.44	8.32	8.21	8.09	7.98
24°C	8.62	8.51	8.40	8.28	8.17	8.06	7.95	7.84
25°C	8.47	8.36	8.25	8.14	8.03	7.92	7.81	7.70
26°C	8.32	8.21	8.10	7.99	7.89	7.78	7.67	7.56
27°C	8.17	8.07	7.96	7.86	7.75	7.64	7.54	7.43